

Your new route to independence

# Rascal Vista & Vista DX Profile:

Introducing the exclusive Rascal Vista and Vista DX mid-size pavement scooters. Benefiting from the quick release feature it dismantles into a few easy-to-lift sections for transporting in your car. The compact & slim-line style offers plenty of legroom yet it is still extremely robust with generous turning manoeuvrability, making it ideal for everyday use.

Each available in two eye catching colours with appealing features including easy access tiller charging & delta handle bar plus the all-round suspension offers added comfort. When taking apart, the split battery pack is also a real advantage. It comes with front & rear LED lights maximising visibility, giving you both confidence & support wherever you go. Available with solid puncture proof or pneumatic tyres - the choice is yours!



### Enjoy day trips without the worry!

The Vista DX's increased 50Ah battery size offers an impressive range of up to 34.5 miles & the luxury Captain's seat makes it an excellent solution for active people who require a multi-purpose scooter & want to stay out for longer with friends & family.

Vista DX shown in Sunset Orange



#### Main Features Vista:

- All-round suspension
- Easy take apart action car transportable
- Dual compound solid puncture proof tyres
- Comfort seat with folding backrest, adjustable arms & integrated swivel function
- · Delta tiller with useful charging port
- · Convenient front basket
- · Class leading weight capacity
- · Off-Board Battery Charging capability

# Main Features Vista DX:

- · All of the above features PLUS
- Captain's seat with folding backrest, headrest, adjustable arms with swivel & slide functions
- Range of up to 34.5 miles with 50A Battery capacity



Travel in comfort & style with the all-round suspension and slide & swivel seating





OPTIONAL

## Options:

- · Rascal Armrest bag
- Rear Saddle bag
- · Metal Cane holder kit



Also available with Pneumatic tyres

#### Fully dismantles into 6 easy to handle parts for storage & transportation




Tyre options & colour choices available on both models

# Vista & Vista DX

Specification:	Vista	Vista DX
Type Class	B/2	$\checkmark$
Overall length	118cm / 46.5"	✓
Overall width	58.5cm / 23"	$\checkmark$
Floor to seat height	57cm - 62cm	$\checkmark$
Seat width	44.5cm / 17.5"	45.7cm / 18"
Seat depth	43.1cm / 17"	43.1cm / 17"
Maximum carrying capacity	150kgs / 23.6 stone	140kgs / 22 stone
Battery voltage and capacity	12V / 35Ah x 2	12V / 50Ah x 2
Maximum speed	6.4km/h / 4mph	$\checkmark$
Ground clearance	7.5cm / 3"	$\checkmark$
Maximum kerb height	6cm / 2.3"	$\checkmark$
Range*	up to 38km / 24 miles	up to 55.5km / 34.5 miles
The maximum safe slope	6°	$\checkmark$
Turn-around width	1.7m / 5'8"	$\checkmark$
Colour	Oxford Blue / Sunset Orange	Oxford Blue / Sunset Orange
Single Heaviest component	25kgs / 55.12 lbs	$\checkmark$

\*Range on full charge and flat ground See full Tech Page on Website & Owner's Manual for details



(front Chassis)





Local Rascal Mobility Retailer:



