

ELECTRIC
MOBILITY®



The Ambassador

 **BRITISH
MADE**
Quality Craftsmanship

Tilt-in-Space Dual Motor




Waterfall Back



Alta Mocha



Alta Charcoal



- Leading global producer of high quality, finished leather
- Natural material
- Quality leather furniture will have a life span of many years

Featuring Variable Angle Lift*


CosiChair




Images are for illustrative purposes only.



The Ambassador

Waterfall Back

The Ambassador Riser Recliner in Leather is the perfect example of quality British craftsmanship combined with timeless, durable and stylish design. The deep waterfall cushions provide exceptional standards of comfort and relaxation. The seating position offers ultimate support even while rising or reclining. Make the right choice and invest in a premium leather chair for the long term.



crestleather


Leather made using small amounts of finishing can outlast fabric upholstery material by between four and five times.

Produced with minimum environmental impact.

Size Guide	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg

* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

Main Features & Benefits:

- Seminar Dual-motor Tilt-in-Space Mechanism with Variable Angle Lift* (VAL* - M & L only. Ask for more details)
- Hybrid pocket sprung seat base
- Pleated wings for a timeless elegance
- Arm covers and Antimacassar included
- Heat & massage available Made to order
- Other fabrics are available on Made to order
- 2 or 3-Seater static sofa available on Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- **Quality Craftsmanship, BRITISH MADE** 
- Leather Care Kits available



Riser Recliners available from **Stock** in the following fabric



Altara Mocha



Altara Charcoal

Riser Recliners available **Made to Order** in the following fabric




Altara Desert




Altara Adriatic



Altara Sage



Altara Smoked Oyster



Altara Mushroom



Altara Chocolate



Altara Hellbraun

