

Vie

Transportable Scooters



'Best Value Boot Scooter'



Affordability for **Everyone!**

RASCAL



Vie Specification:

Type Class	B / 2
Overall length	108.3cm / 42.6"
Overall width	47cm / 18.6"
Floor to seat height	39-43cm (min-max)
Seat width	41cm / 16.1"
Seat depth	37cm / 14.6"
Maximum carrying capacity	121kgs / 19 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	4.5cm / 1.8"
Maximum kerb height	5cm / 2"
Range* (Max kg user weight)	up to 14.2km / 8.8 miles
The maximum safe slope	6°
Turn-around width	1.25m / 4'1"
Colour	Blue / Red

Single Heaviest component 14.8kgs / 32.5lbs

*Range on full charge and flat ground. Tested at maximum carrying capacity.
See full Tech Page on Website & Owner's Manual for details

Rascal Vie Profile:

Are you ready for the Rascal Vie? Here comes affordability for everyone. Economically priced, with a range upto 8.8 miles makes it suitable for your daily run about. This 4mph scooter is comfortable to use, yet rugged & nimble with a tight turning circle. Transportable with the front section weighing just 14.8kg, giving you all round versatility & daily enjoyment when you are outdoors with friends & family.

The Rascal Vie is a very practical scooter with an easy & robust take apart mechanism which allows you to store or drive away with it in the boot of your car, perfect for everyday. Experience this class-leading boot scooter for yourself.

Local Rascal Mobility Retailer:

Main Features:

- **A quality product, economically priced**
- Manoeuvrable with tight turning circle
- Car transportable
- Delta tiller & convenient charging point
- Comfort swivel seat with white stitching
- Solid foam filled puncture proof tyres
- Spacious shopping basket
- Rear safety reflectors

OPTIONAL

Optional:

- Rascal arm rest bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v3 July 2025