

Veo Sport SR

Transportable scooters



Your new route to independence

RASCAL

Veo Sport SR

Veo Sport LiFe Model with
Lithium Ferro-Phosphate
Battery also available

Veo Sport SR Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	39-44cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	130kgs / 20.5 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.4"
Maximum kerb height	5cm / 2"
Range* (Max kg user weight)	up to 25.1km / 15.6 miles
Range* (90kg user weight)	up to km 29.2/ 18.2 miles
The maximum safe slope	6°
Turn-around width	1.40m / 4'6"
Colour	Neon Blue / Mulberry
Single Heaviest component	20.2kgs / 44.5lbs

*Range on full charge and flat ground.

See full Tech Page on Website & Owner's Manual for details

Rascal Veo Sport SR Profile:

It's all about the drive. The Rascal Veo Sport SR is fully loaded with features in 2 striking colours. A modern, sporty design with class leading all-round suspension. With 22Ah batteries, a superb capacity for a scooter of its size, affording a range of up to 15.6 miles. A padded delta tiller boasting a convenient charging point and luxury seat, detailed with unique white accent stitching, the Veo Sport SR is a real pleasure to ride.

The Veo Sport SR's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension plus LED front light, this is the travel scooter with the edge.

Local Rascal Mobility Retailer:

Main Features:

- Car transportable
- All-round suspension
- Stylish sports seat
- **Tiller charging point**
- **LED front light**
- **Front bumper**
- Dual compound puncture proof tyres
- Comfortable delta tiller
- Rear safety reflectors



OPTIONAL



Optional:

- Rascal arm rest bag



The LED front light & bumper together with class leading all-round suspension unites design, functionality & safety.



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v5 August 2025

