

ELECTRIC
MOBILITY®



 **BRITISH
MADE**
Quality Craftsmanship

The Jubilee

Tilt-in-Space Dual Motor



Lateral & Waterfall Back



Featuring Tilt-in-Space

and Variable Angle Lift*


CosiChair





The Jubilee

Lateral & Waterfall Back


The Jubilee is our UK Premium Tilt-in-Space (TIS) Dual-motor Cosi Chair. Designed to maximise support throughout motion with a range of stylish back designs, upholstery fabrics and offering a high leg elevation to provide exceptional postural seating positions to suit most individual requirements. The Jubilee is available in Lateral and Waterfall Back from stock, with other fabrics available Made to Order.



Size Guide	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg

Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to Order

Main Features & Benefits:

- Dual-motor Tilt-in-Space mechanism with Variable Angle Lift (VAL* - Medium and Large only. Ask for more details)
- Hybrid pocket sprung seat base
- Arm covers and Antimacassar included
- 5-button handset
- Choice of placement with configurable back cushions
- Heat and massage available Made to Order
- Drop in seat cushions for pressure care issues available Made to Order. (Ask for more details.)
- 2 and 3-Seater static sofas available on Made to Order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, BRITISH MADE 

Exclusively available from your local Cosi Chair Specialist




Café



Duck Egg



Forest Green



Marine



Slate



Stone



Mink



Oatmeal



Sky



Plum



Cocoa



Zinc

