

Vista LiTEWAY

Lithium powered scooters

**LiFetime Battery
Guarantee***



**RASCAL
Life**



The Lithium Ferro-Phosphate Battery
is built to last and is one of the safest
Lithium Batteries in the Mobility sector

RASCAL



Vista LiTEWAY

LiFetime Guarantee
on Lithium Ferro-Phosphate Battery†

† T&Cs apply

Vista LiTEWAY Specification

Type Class	B / 2
Overall length	118cm / 46.5"
Overall width	58.5cm / 23"
Floor to seat height	55cm - 60cm
Seat width	48cm / 18.9"
Seat depth	42cm / 16.5"
Maximum carrying capacity	140kgs / 22 stone
Battery voltage & capacity	25.6V / 27Ah
Maximum speed	6.4km/h / 4mph
Ground clearance	6.8cm / 2.68"
Maximum kerb height	5cm / 2"
Range* (Max kg user weight)	up to 26km / 16.3 miles
Range* (90kg user weight)	up to 36km / 22.5 miles
The maximum safe slope	6°
Turn-around width	1.58m / 5'2"
Colour	Aston Green / Oxford Blue
Single Heaviest component (front chassis)	14.7kgs / 32.4 lbs

*Range on full charge and flat ground.

See full Tech Page on Website & Owner's Manual for details.

Main Features:

**RASCAL
Life**

- 27Ah **Lithium Ferro-Phosphate Battery**
- All-round suspension
- Easy take apart action - transportable
- Generous padded seat with folding backrest, adjustable arms and integrated swivel function
- Sleek, lightweight aluminium frame
- Delta tiller with useful charging port and increased height for upper body comfort
- Off-board Battery Charging capacity

The choice is yours:

Available with solid
puncture proof or
pneumatic
tyres!



Rascal Vista LiTEWAY Profile:

Introducing the Rascal Vista LiTEWAY with a LiFetime Battery Guarantee. This mid-size aluminium scooter benefits from the quick release feature and it dismantles into a few easy-to-lift sections for transporting in your car. The compact and slim-line style offers a spacious footwell and plenty of legroom, yet it is still extremely robust with generous turning manoeuvrability, making it ideal for everyday use.

Available in two eye catching colours with appealing features including easy access tiller charging and delta handle bar plus the all-round suspension offers added comfort. It comes with front and rear LED lights maximising visibility, giving you both confidence and support wherever you go.

BATTERY
RASCAL
Life

27Ah battery, weighing just 7.2kg:

Relax and enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO4) Battery pack is the safest, **greenest** with an increased efficiency and performance to stand up to everyday life.



† LiFetime Battery Guarantee T&Cs apply



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v3 June 2026

 electricmobility.co.uk