

# Pioneer

Performance scooters



Your new route to independence

**RASCAL**

# Pioneer



## Pioneer Specification:

Type Class	C / 3
Overall length	147cm / 58"
Overall width	64.5cm / 25"
Floor to seat height	65-68cm (min-max)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity	12V / 75Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range* (Max user weight)	up to 43.7km / 27.3 miles
Range* (90kg user weight)	up to 56km / 35 miles
The maximum safe slope	10°
Turn-around width	2.02m / 6'7"
Colour	Matt Blue / Gloss Red

\*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

## Rascal Pioneer Profile:

The Rascal Pioneer, a proven performer with fantastic reliability and durability record has had a style refresh. With full coil suspension and a 35 mile range (90kg user weight) its powerful motor makes it an ideal scooter for the more challenging environments. It features a sliding swivel seat with headrest and padded arms, showerproof switches, a large basket and a practical handy storage compartment for charging your phone on the move. Choose the Rascal Pioneer, one of the most popular scooters on the market.

- **Tried & Trusted**
- **Simple & Easy to Use**

## Main Features:


- LED front light
- Front and rear coil suspension
- Practical front basket
- Gas strut tiller adjustment with useful battery charging port
- Large shock absorbing tyres for maximum comfort on the move
- Sliding swivel seat with armrests and headrest for maximum comfort

## OPTIONAL ACCESSORY



## Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- An exclusive Rascal basket bag
- Puncture proof infill tyres



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E. - v12 February 2026