

Technical Specification - Rascal ROLL8 21 & 23 Rollator

Rollator Type	ISO 11199-2: 2021	Outdoor & Indoor
Max. Dimensions fully assembled	length x width x height mm (in)	720(28.3) x 605(23.8) x 828(32.6) - ROLL8 21 720(28.3) x 605(23.8) x 914(36.0) - ROLL8 23
Min. Dimensions- folded & handles retracted	length x width x height mm (in)	720(28.3) x 230(9.0) x 682(26.8) - ROLL8 21 720(28.3) x 230(9.0) x 737(29.0) - ROLL8 23
Maximum carrying capacity	Kg (lbs) ((stone))	130 (286) ((20.5))
Mass (weight)	Kg (lbs)	7.2 (15.8) - ROLL8 21 7.3 (16.0) - ROLL8 23
Seat Width	mm (in)	390 (15.4)
Seat Depth	mm (in)	205 (8.1)
Seat Height	mm (in)	550 (21) - ROLL8 21 580 (23) - ROLL8 23
Rear Wheel Size	mm (in)	230 x 42 (9 x 1.6)
Front Castor Wheel Size	mm (in)	230 x 42 (9 x 1.6)
Type of wheels/ tyres	Type	PU
Storage Bag Maximum Load	Kg (lbs)	5 (11)

NOTE: This Rollator meets the relevant requirements of ISO 11199-2: 2021

Due to a policy of continual improvement, Electric Mobility Euro Ltd. reserves the right to change specifications without prior notice.