

Rascal Roll8

Rollator



Your new route to independence

RASCAL



Rascal Roll8



Roll8 Specification:

Handlebar height (min & max)	77cm (30.3") / 90cm (35.4") (Roll8 21) / 77cm (30.3") / 93cm (36.6") (Roll8 23)
Max. dimensions fully assembled (L x W x H)	72cm (28.3") x 60.5cm (23.8") x 82.8cm (32.6") (Roll8 21) / 72cm (28.3") x 60.5cm (23.8") x 91.4cm (36") (Roll8 23)
Min. dimensions folded & handles retracted (L x W x H)	72cm (28.3") x 23cm (9") x 68.2cm (26.8") (Roll8 21) / 72cm (28.3") x 23cm (9") x 73.7cm (29") (Roll8 23)
Seat height	55cm / 21" (Roll8 21) / 58cm / 23" (Roll8 23)
Seat width	39cm / 15.4"
Seat depth	20.5cm / 8.1"
Maximum carrying capacity	130kgs / 20.5 stone
Wheel size	23cm x 4.2cm / 9" x 1.6"
Colour	Black
Storage bag maximum load	5kgs / 11lbs
Mass (Weight)	7.2kgs / 15.8lbs (Roll8 21) / 7.3kgs / 16lbs (Roll8 23)

Rascal Roll8:

Lightweight. Stylish. Exceptionally stable. Crafted from high-performance carbon, this rollator delivers not only strength and durability but also a sleek contemporary design to help you walk taller. With 'easy on the hands' ergonomic handlebars for a stable grip and 9" soft roll PU wheels, every step feels smooth and steady. Use the clever kerb assist to go up pavements without worry and enjoy the reassurance of the Rascal Roll8 rollator. 2 Year Frame warranty.



Main Features:

- **Lightweight Design**
At just 7.2kg, this rollator is incredibly easy to lift, transport and store, ideal for both users and carers.
- **Two Size Options for the Perfect Fit**
Available in 21" and 23" seat heights, ensuring improved posture and personalised comfort.
- **Superior Stability and Control**
The 9" wheels on all four corners provide enhanced balance and smooth manoeuvrability across indoor and outdoor surfaces.
- **Convenient Storage**
Detachable front bag with 5kg capacity keeps essentials secure and easily accessible.
- **Comfortable and Supportive Seating**
Built-in seat supports up to 20.5 stone / 130kg, offering a safe and comfortable resting point whenever needed.

Height adjustable ergonomic handlebars



Regain confidence in every step, whether it's a trip to the shops or a walk outdoors, this rollator is designed to support independence, comfort and peace of mind.

[A trusted Rascal.](#)