

Rascal C60 Carbon

Compact & folding powerchairs



**LiFetime Battery
Guarantee***

The Lithium Ferro-Phosphate Battery is built to last and is one of the safest Lithium Batteries in the Mobility sector



RASCAL



Rascal C60 Carbon

LiFetime Guarantee
on Lithium Ferro-Phosphate Battery*

† T&Cs apply

C60 Carbon Specification:

Type Class	B / 2
Overall length	91cm / 35.8"
Overall width	64cm / 25.2"
Floor to seat height	50cm
Seat width	41cm / 16"
Seat depth	43cm / 17"
Maximum carrying capacity	120kgs / 18.9 stone
Battery voltage and capacity	24V / 10Ah x 1
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	4cm / 1.6"
Range* (Max user weight)	up to 13km / 8.1 miles
Range* (100kg user weight)	up to 15km / 9.4 miles
The maximum safe slope	6°
Turn-around width	1.2m / 3'9"
Colour	Black
Single Heaviest component	14.3kg / 31.5lbs

*Range on full charge and flat ground

*See full Tech Page on Website & Owner's Manual for details

Rascal C60 Carbon:

The Rascal C60 carbon powerchair combines cutting-edge technology with user-friendly design. Featuring all-round suspension and a high back seat with Memory Foam cushion for optimum support and comfort for everyday use. A brushless motor offers a smooth, quiet ride delivering responsive acceleration. The quick release battery is practical for travel. The C60 Carbon Powerchair designed for independence and performance.



Main Features:

- **Ultimate portability, worry-free travel**
The Rascal C60 Carbon achieves the perfect balance between lightweight and durability. Its sophisticated foldable design makes storage easy, allowing you to travel anywhere.
- **Long-lasting battery life**
The Rascal C60 Carbon is equipped with a **high-performance Lithium Ferro-Phosphate battery**, with Rascal LiFetime Battery Guarantee for peace of mind.
- **Smart control, easy control**
The Rascal C60 Carbon's joystick control system makes operation simple. Ergonomically designed to prevent fatigue even after extended use.
- **Safe and secure travel**
Equipped with multiple safety features, including anti-roll wheels and an intelligent braking system, ensures safe travel at all times.



OPTIONAL ACCESSORY

- Bluetooth Remote Controller



Air travel & CAA Lithium Battery Legislation



Airline Compliant 10Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel and **can**

be taken on a plane.
(Subject to Airline Regulation and Restrictions)

Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v1 June 2026