

Vecta Cruise



RASCAL 

Performance Scooter

Vecta Cruise

The Vecta Cruise is a 8mph performance scooter that combines design with technology to make a scooter that looks good and performs even better. The Vecta Cruise is robust and sturdy on the road, agile enough for roads and pavements.

Main Features:

- Range up to 27.3 miles*
- Up to 25 stone carrying capacity
- “Max Grip” transmission, and featuring a limited slip differential gearbox
- 600 Watt high power motor
- Premium Captain’s seat with elegant diamond-quilted stitching, lumbar support and flip up armrests for easy access
- Super bright LED headlights with smoke tinted lenses
- Adjustable tiller and seat for optimum ergonomic positioning
- 52 inches long
- Pneumatic tyres
- **Tevere Green**





Optional upgrades:

- Solid Canopy
- Solid Canopy with Textile Sides
- Rear Top Box



18" wide highly-padded seat with lumbar support and flip up armrests



Adjustable seating

25-stone carrying capacity



Pneumatic tyres



High torque 600 Watt motor and "Max Grip" transmission

52 inches long

Mirrored running light

Range up to 27.3 miles*

Adjustable tiller

Storage pocket and a dedicated USB port



Smoke tinted lenses

Sleek and streamlined

Vectability

Vecta Cruise

Vecta Cruise Specification

Type Class	C / 3
Overall length	132cm / 52"
Overall width	60cm / 23.6"
Floor to seat height	62 - 67cm (min-max)
Maximum carrying capacity	160kgs / 25 stone
Battery voltage and capacity	12V / 50Ah x2
Maximum speed	12.87km/h / 8mph
Ground clearance	10cm / 4"
Maximum kerb height	10cm / 4"
Maximum safe descendable kerb height	10cm / 4"
Range*	up to 44km / 27.3 miles
The maximum safe slope	10°
Turn-around width	1.7m / 5'5"
Turning diameter	3.47m / 11'38"
Colour	Tevere Green
Lighting	LED front and rear
Wheel diameter	30.5cm / 12"
Type of tyres	Pneumatic

*Range on full charge and flat ground

*See full Tech Page on Website & Owner's Manual for details

Local Rascal Mobility Retailer:



Max Grip-Limited Slip Differential Gearbox



600 Watt High Power Motor



Suspension



USB Port



LED



Vecta Cruise with and without Canopy meets the requirements as a Class 3 Scooter with necessary user equipment as defined with the Invalid Carriage Act (Amended 2015)



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O. E- v1 July 2026