

Vortex Cruise



RASCAL 

Performance Scooter

Vortex Cruise

The Vortex Cruise is a unique and dynamic 8mph scooter ready to take your mobility and lifestyle to the next level. Combining design with the latest technology, the Vortex Cruise is built to perform with a host of user-friendly features.

Main Features:

- Range up to 25.5 miles*
- Up to 28 stone carrying capacity
- Premium 20" wide Captain's seat with elegant diamond-quilted stitching, lumbar support and flip up armrests
- **Racing Green**





2 rear view mirrors

Mirrored running light

Spacious legroom

Adjustable tiller and seat
for optimum ergonomic
positioning

900 Watt high torque motor
and Max Grip limited slip
differential gear box

Pneumatic tyres



Super bright LED lights
with smoke tinted lenses

Damped suspension

Waterproof USB
charging point

Optional upgrades:

- Solid Canopy
- Solid Canopy with Textile Sides
- Rear Top Box



Travel in first-class comfort

Vortex Cruise





Vortex Cruise Specification

Type Class	C / 3
Overall length	146cm / 57.5"
Overall width	65cm / 25.6"
Floor to seat height	62-67cm (min-max)
Maximum carrying capacity	180kgs / 28 stone
Battery voltage and capacity	12V / 75Ah x2
Maximum speed	12.87km/h / 8mph
Ground clearance	14cm / 5.5"
Maximum kerb height	10cm / 4"
Maximum safe descendable kerb height	10cm / 4"
Range* 75Ah battery	up to 41km / 25.5 miles*
Range* 100Ah battery	up to 54km / 33.5 miles
The maximum safe slope	10°
Turn-around width	2.2m / 7'2"
Turning diameter	4m / 13'01"
Colour	Racing Green
Lighting	LED front and rear
Wheel diameter	34cm / 13.4"
Type of tyres	Pneumatic

*Range on full charge and flat ground
See full Tech Page on Website & Owner's Manual for details

Local Rascal Mobility Retailer:

Vortex Cruise 75Ah with and without Canopy meets the requirements as a Class 3 Scooter with necessary user equipment as defined with the Invalid Carriage Act (Amended 2015)



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O. E- v1 July 2026

electricmobility.co.uk