

Juno

Transportable Scooters



RASCAL

Affordability for **Everyone!**

Juno

Juno Specification:

Type Class	B / 2
Overall length	103.5cm / 40.7"
Overall width	49.5cm / 19.5"
Floor to seat height	45-49cm (min-max)
Seat width	43cm / 17"
Seat depth	35.5cm / 14"
Folded Dimensions LxWxH	103.5cm / 49.5cm / 71cm (40.7" / 19.5" / 27.9")
Maximum carrying capacity	136kgs / 21.4 stone
Battery voltage and capacity	12V / 14Ah x2
Maximum speed	6.4km/h / 4mph
Ground clearance	4.5cm / 1.8"
Maximum kerb height	5cm / 2"
Range* (Max kg user weight)	up to 13.5km / 8.4 miles
Range* (90kg user weight)	up to 16.5km / 10.3 miles
The maximum safe slope	6°
Turn-around width	1.25m / 4'1"
Colour	Graphite Grey / Steel Blue / Diamond White

*Range on full charge and flat ground.

See full Tech Page on Website & Owner's Manual for details

Rascal Juno Profile:

Discover the Rascal Juno. An economically priced and versatile 4mph transportable scooter. With a range of up to 10.3 miles (90kg user weight) makes this ideal for your daily run about. Comfortable with a soft padded seat and suspension, its designed to encourage good posture with the increased legroom. Nimble with a tight turning circle, the Juno is engineered for stability for everyday errands. A premium finish with practical features including a USB port and tiller lock for easy lifting and storage when dismantled.

Main Features:

- **A reliable and quality Rascal**
- Manoeuvrable with tight turning circle
- Solid foam filled puncture proof tyres
- Increased legroom, suitable for a range of heights from 4'11" – 6'4"
- Spacious shopping basket
- Suspension
- Tiller lock
- Practical USB port



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v1 July 2026