

ALLTRACK N21

NARROW 21" BASE MWD

EXTRA-NARROW BY DESIGN. POWERFUL BY NATURE.

When space is limited but performance is non-negotiable, the Alltrack Narrow N21 delivers. With a slim 21" base width, this model was designed to navigate tight corners, narrow doorways and compact living environments without compromising on comfort, stability or functionality. The N21 offers our patented Dual-Action Suspension for smooth travel indoors and out. Whether you're moving through the home, the office or around town, this award-winning powerchair gives you the confidence to go wherever life takes you.

Highlights

- The N21 features a 21" base width with 14" drive wheels
- True MWD design provides a compact 20" turning radius up to 1" smaller than other narrow bases
- Seat sizes from 14" x 14" up to 22" x 22"
- Seat-to-floor height as low as 15.5"



21"

Shown in Candy Red

RASCAL
REHAB


ELECTRIC
MOBILITY

Exclusive Distributor

Technical Specifications

Configuration	Mid-Wheel Drive with Narrow Base (Alltrack N21)
Suspension system	Patented Dual-Action Suspension: Fully independent 6-wheel interactive suspension with vertical shock absorption
Range up to	30 miles / 48.2 km/h
Ground clearance¹	3.5" with 14" drive wheels, 2.875" with 12.5" drive wheels
Base width²	N21 base width 21" with 55Ah batteries
Base length²	35.125"
Turning radius	20"
Weight capacity	21.4 stone (136 kg)
4-Pole motor options	6.5 mph (10 km/h)
Transport option	4-point tie-down attachment points compliant with ISO 7176-19
Battery types	Two 12 V, deep cycle Gel cell Group 22NF (50 Ah)
Battery charger	8 A, off board
Electronics	VR2 (non-expandable 90 A) and R-net (expandable 90 A or 120 A)
Specialty controls	Compatible with Mo-Vis, ASL, Switch-It, Munevo and Stealth as well as other current systems
Seat width range	Adjustable from 14" to 22" wide
Seat to floor height range (Tilt <u>OR</u> Seat Elevate)	15.5" to 19.5"
Seat to floor height range (Tilt <u>AND</u> Seat Elevate)	17.5", 18.5" or 19.5"
Drive wheel sizes	12.5" x 2.25" or 14" x 3"
Total weight	199.2lb (90.2 kg)
Castor wheel size	6"

1. Subject to variances depending on configuration and manufacturing tolerances.
2. Dimensions listed are for power base only.



Exclusive Distributor

Electric Mobility contact details:

Tel: 01460 258120

Email: rehab@electricmobility.co.uk

Amylior - We reserve the right to make technical changes and cosmetic changes to our products
Images are for illustrative purposes only
E. and O.E- v2 September 2026

 electricmobility.co.uk



Electric Mobility
Canal Way
Ilminster
Somerset
TA19 9DL

