

ALLTRACK^{R3} HYBRID

HYBRID-WHEEL DRIVE (HWD) CUSTOM OFF-ROAD POWERCHAIR

**NEXT-LEVEL OUTDOOR PERFORMANCE.
NO COMPROMISE INDOORS.**

No Compromise! – The award-winning Alltrack R3 Outdoor is engineered for users who require indoors manoeuvrability but need a chair to handle challenging rural and/or urban environments. With its one of a kind HWD footprint, extra-wide tyres, larger castors and advanced suspension, it delivers unstoppable traction and control over uneven, unpredictable terrain. Whether it's dirt paths, gravel tracks, or snowy trails, this powerhouse keeps you moving forward with confidence. Despite its off-road capabilities, the R3 Outdoor still delivers the compact agility needed indoors – making it the ultimate off-road companion for everyday adventurers.

Highlights

- Outdoor-ready Hybrid-Wheel Drive Technology for off-road performance and indoor manoeuvrability
- Award-winning, rugged design tested by real users in real environments
- Powerful 4-pole motors
- 9" front and rear castors offer excellent grip and enhanced ride comfort
- Standard 14"x 4" drive wheels for added comfort and grip on loose, soft terrain
- Flared front-end for enhanced castor clearance and added stability
- Optional Hydrographic finishes to match your style



Shown in Woodland Camo Hydrographic

mobility ★
product award 2019

RASCAL
REHAB




**ELECTRIC[®]
MOBILITY**

Exclusive Distributor

Technical Specifications

Configuration	Hybrid-wheel drive – Standard (Alltrack R3 Hybrid with Off-Road Powerchair)
Suspension system	Patented Dual-Action Suspension: Fully independent 6-wheel interactive suspension with vertical shock absorption
Range up to	26 miles / 41.8 km/h (55Ah batteries) 37 miles / 59.5 km/h (70Ah batteries)
Ground clearance¹	3.5" with 14" drive wheels
Base width²	23.375" with 55Ah batteries or 25.375" with 70Ah batteries
Base length²	40.9"
Turning radius	21.75"
Weight capacity	21.4 stone (136 kg)
4-Pole motor options	6.5 mph (10 km/h) / 8.5 mph (13.5 km/h) ³
Transport option	4-point tie-down attachment points compliant with ISO 7176-19
Battery types	Two 12 V, deep cycle Gel cell Group 22NF (55 Ah) or Group 24 (70 Ah)
Battery charger	8 A, off board
Electronics	R-net (expandable 120 A)
Assistive driving package	Available with Smart-Track technology
Specialty controls	Compatible with Mo-Vis, ASL, Switch-It, Munevo and Stealth as well as other current systems
Seat width range	Adjustable from 14" to 26" wide
Seat to floor height range (Tilt <u>OR</u> Seat Elevate)	16.75" to 19.75"
Seat to floor height range (Tilt <u>AND</u> Seat Elevate)	18.75" or 19.75"
Total weights	55Ah version: 199.2 lb (90.2 kg) 70Ah version: 240 lb (108 kg)
Drive wheel size	14" x 4"
Castor wheel size	9"

1. Subject to variances depending on configuration and manufacturing tolerances.
2. Dimensions listed are for power base only.
3. Important UK usage notice: This powerchair has a maximum speed capability of 8.5 mph. Users must not exceed 8 mph on public roads or 4 mph on pavements and in pedestrian areas. The user is responsible for ensuring that the powerchair is registered, insured where required, and operated in accordance with applicable UK legislation



Exclusive Distributor

Electric Mobility contact details:

Tel: 01460 258120

Email: rehab@electricmobility.co.uk

Amylior - We reserve the right to make technical changes and cosmetic changes to our products
Images are for illustrative purposes only
E. and O.E- v2 September 2026

 electricmobility.co.uk



Electric Mobility
Canal Way
Ilminster
Somerset
TA19 9DL

