



Redefining Air Cushion Technology



Made to Measure
available in 2", 3" and 4" cell heights



Exclusive Distributor

Performa CG Air Cushion



CG Air Cushion Profile:

Designed for superior comfort, stability and skin protection, these advanced cushions use cutting edge materials to provide an unparalleled combination of support and pressure reduction aimed at delivering better overall patient outcomes.

Through advancements in materials and manufacturing techniques, Performa CG Air is a true reinvention of the air cushion.

Designed to tackle industry-wide challenges - including hygiene, stability, durability, maintenance, and positioning - CG Air delivers unmatched performance without compromising pressure management.

With its enhanced benefits, Performa CG Air opens new possibilities for users seeking the best in pressure management therapy.

Main Features:

Medical Grade PVC Cell Design

- Soft to the touch material offers first class stability and durability
- Resistant to discolouration and odour



Cushion Configurations

- Single or dual compartment configurations
- Custom dimensions and cell configurations, asymmetrical seat depths available

Breathable Cover Choices and Advancements

- Covers are made of spacer knits in single layer (standard) and 3D double layers (optional), which reduce shear, enhance comfort, and promote air flow.
- GripZone Patent Pending anti-slip technology (Optional) - *Patent Pending
- Prevents shear forces by placing a thin anti-slip polymer under the femurs and trochanters while protecting sensitive areas.
- No reduction in fabric's stretch or breathability
- Stitch-free design



Standard cover



Cover with GripZone cover

Superior Pressure Management

- CG Air's innovative air channels provide exceptional pressure reduction by providing a smoother, controlled airflow beneath the user, helping to stimulate blood circulation and reduce pressure points.

Enhanced Postural Stability and Positioning

Made of soft to the touch low stretch, medical grade PVC, CG Air retains its shape even if overinflated, providing a more stable seating surface. Unlike traditional air

cushions, the unique cell material actively creates an ischial well and helps to reduce pelvic migration and posterior pelvic tilt, creating an unparalleled combination of pressure reduction and stable positioning.

Better Air Retention and Hygiene

CG Air retains proper inflation and is highly resistant to leaks due to its patented High-Frequency welded design and use of stainless-steel valves instead of plastic.

The medical grade PVC is also easy to maintain daily with soap and chlorine solutions and is even machine washable and autoclave safe.

Dermatologically friendly

- Crafted from fire/flame-retardant (CAL117-2013), non-allergenic, and anti-fungal medical grade PVC, CG Air is ideal for users with allergies or specialised care needs.

Optional extras:

Anatomic Contour Seat Insert -
increases immersion and lateral stability



Transfer Trellis -
supports transfers while protecting vulnerable front cells



Convenient Seat Pouch -
Perfect for a mobile phone and everyday essentials with easy magnetic flap



Exclusive Distributor



Electric Mobility
Canal Way
Ilminster
Somerset
TA19 9DL

Amylior - We reserve the right to make technical changes and cosmetic changes to our products
Images are for illustrative purposes only
E. and O.E- v2 September 2026

Tel: 01460 258120

Email: rehab@electricmobility.co.uk

 electricmobility.co.uk