

# Ultralite 480

Transportable scooters



Your new route to independence

**RASCAL** 

# Ultralite 480

Voted "Best Buy" in the Mobility Scooter Reviews 2014 by **Which?** Magazine

## UL480 Specification:

Type Class	A / 2
Overall length	102cm / 40"
Overall width	51cm / 20"
Floor to seat height	45 - 55cm
Seat width	43.1cm / 17"
Seat depth	33cm / 13"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	4cm / 1.5"
Range* - standard battery pack	up to 13km / 8 miles
The maximum safe slope	6°
Turn-around width	1.62m / 5'4"
Colour	Black / Red / Teal
Single Heaviest component	13kgs / 29lbs

\*Range on full charge and flat ground

## Rascal Ultralite 480 (UL480) Profile:

The ultimate transportable scooter, with its compact & lightweight design, the Rascal UL480 ensures you stay independent at all times. This scooter has a proven reputation for outstanding design & reliability. Easy to transport, the UL480 can be dismantled in seconds; no cables to worry about, just quick release locks, convenient lifting handles and lightweight manageable components. Easy to stow it can be stored vertically in your home if space is an issue, & fits easily into most car boots when taken apart.

Local Rascal Mobility Retailer:

## Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Swivel seat

**OPTIONAL  
ACCESSORY**



## Optional Extras:

- Optional Rascal arm rest bag
- Cane & Crutch Holder



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E - v8 February 2019