Rascal Luxe

Lithium powered scooters



LiFetime Battery Guarantee*



RASCAL LiFe

The Lithium Ferro-Phosphate Battery is built to last and is one of the safest Lithium Batteries in the Mobility sector

40 Years ANNIVERSARY CELEBRATION

RASCAL EXCLUSIVE

RASCAL EXCLUSIVE

LiFetime Battery Guarantee†


† T&Cs apply

Rascal Luxe Profile:

The **Rascal Luxe** 4mph **Lithium Ferro-Phosphate** powered transportable scooter is sleek and durable. With a lightweight aluminium frame, designed for ergonomic “All day” comfort offering increased legroom that accommodates a range of heights from 4’11” - 6’4” plus.

This premium scooter features a supremely comfortable deluxe seat, encouraging **100°** posture perfect alignment of the feet, legs and back position for a supportive journey. The space around the tiller, offers significant leg and knee room that promotes even weight distribution and enhances the stability alongside the adjustable all-round suspension.





The **Rascal Luxe** benefits from a quality design that showcases styling, an impressive ride up to 21.4 stone user weight and reliable range up to 21.7 miles (90kg user weight). Luxuriously developed comfort and manoeuvrability for every journey with this class leading scooter!

Main Features:



Lithium Ferro-Phosphate Battery

Where exceptional design meets safety...





“See and be seen” with the **360° LED lighting,** highly visible from any angle



Useful **tiller lock** helps with lifting, **weighs just 10.8kg** and super portable for easy storage





Adjustable **all-round suspension** Designed to absorb uneven pavements, reduce impact and user fatigue

IPX5ˣ waterproof rating. ˣ Ingress Protection Rating

Practical **USB port**





The **Rascal Luxe** with high quality **Lithium Ferro-Phosphate** (LiFePO4) Battery. Relax and enjoy the ride in the knowledge that the Lithium Ferro-Phosphate (LiFePO4) Battery pack has an increased efficiency and performance to stand up to everyday life. Plus the **LiFetime Battery Guarantee**ǂ for peace of mind!

Scan to visit website

Truly enhance the experience and freedom to travel

Explorers’ Choice - *Have the Best of Both Worlds*

Airline Compliant

Purchase the 11.6Ah battery option & 2A charger^

Battery weighs just 4.1kg / 9lbs

Range up to 10.3 miles / 16.5km (90kg user weight)*

Air travel & CAA Lithium Battery Legislation

Airline Compliant 11.6Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel and can be taken on a plane. (Subject to Airline Regulation and Restrictions)

ǂ LiFetime Battery Guarantee T&Cs apply

^ Speed charge & save time - charge between 10-90% in less than 2 hours with a 5A charger

Rascal Luxe

LiFetime Guarantee
on Lithium Ferro-Phosphate Battery‡

Specification:

	Luxe 27Ah with Rapid 5A Charger	Luxe 11.6Ah with 2A Charger
Type Class	B / 2	✓
Overall length	103.5cm / 40.7”	✓
Overall width	49.5cm / 19.5”	✓
Overall height	94cm / 37”	✓
Floor to seat height	46-50cm (min-max)	✓
Seat width	44cm / 17.3”	✓
Seat depth	41cm / 16.1”	✓
Folded Dimensions LxWxH	103.5cm / 49.5cm / 71cm (40.7” / 19.5” / 27.9”)	✓
Maximum carrying capacity	136kgs / 21.4 stone	✓
Battery voltage and capacity	25.6V / 27Ah	25.6V / 5.8Ah x 2
Maximum speed	6.4km/h / 4mph	✓
Ground clearance	4.5cm / 1.8”	✓
Maximum kerb height	5cm / 2”	✓
Range* (Max kg user weight)	up to 25.8km / 16 miles	up to 12.2km / 7.6 miles
Range* (90kg user weight)	up to 34.9km / 21.7 miles	up to 16.5km / 10.3 miles
The maximum safe slope	6°	✓
Turn-around width	1.25m / 4’1”	✓
Colour	Champagne / Jade Green / Royal Blue	
Single Heaviest component (front section)	10.8kgs / 23.8lbs	✓
Product weight (Excluding battery) (Including battery)	32kgs / 70.5 lbs 38.7kgs / 85.3 lbs	✓ 36.1kgs / 79.6 lbs

*Range on full charge and flat ground
*See full Tech Page on Website & Owner’s Manual for details

Practical Features:









‡ LiFetime Battery Guarantee T&Cs apply

Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O. E- v6 June 2026