The Haselbury

Dual-Motor Tilt-in-Space

CosiChair By ELECTRIC MOBILITY®









Contour, Lateral and Waterfall Back

The Haselbury

Made to measure

Dual-Motor Tilt-in-Space


The Haselbury Dual-Motor Tilt-in-Space Riser Recliner has been designed with premium, contemporary materials and British quality craftsmanship, giving you the flexibility to choose the size, back style and fabric colour that best suit you and your home. Delivering exceptional functionality, quality and performance for everyday living. With independent dual motors, the backrest and legrest can be adjusted separately, allowing you to find your preferred position that feels right for you. The high leg lift and Tilt-in-Space function works together to provide additional support and comfort by maintaining your seating position as the chair reclines. Where space is limited, this chair features a compact movement with the armrests moving as one in a smooth, gentle action for complete relaxation.
Built to a high standard and backed by a 5 Year Guarantee, the Haselbury is a riser recliner you can enjoy with confidence for years to come.

Size Guide	Small	Medium	Large
Seat height	18½”	20”	21½”
Seat width	19”	21”	22”
Seat depth	18½”	20”	21½”
Arm height	6”	6½”	7”
Back height	27½”	29”	31”
Hip to heel	37½	39”	41”
Overall width	32”	33”	34”
Overall height	43’’	46”	48½”
Weight limit	25 stone	25 stone	25 stone



Knuckles in Light Oak, Teak and Mahogany





Ash



Clay



Espresso



Sage



Easy to use handset

Main Features & Benefits:

- Dual-motor Tilt-in-Space Mechanism
- Useful USB A and C charging ports
- Hybrid pocket sprung seat base
- Pleated wings for a timeless elegance
- Antimacassar and arm covers included
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, BRITISH MADE
- Heat and massage available Made to Order
- Medium Low Base available Made to Order
- Drop in seat cushions for pressure care issues available Made to Order. (Ask for more details.)
- Knuckles available Made to Order
- Other fabrics are available on Made to Order
- 2 and 3-Seater static sofas with fixed back available on Made to Order

CosiChair
electricmobility.co.uk

Images are for illustrative purposes only. Colours may vary slightly from the images shown.
Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice. Details are correct at time of going to print.

E. and O.E - v1 November 2026